

BOLT SYSTEMS, INC.

1700 SILVER STAR ROAD
ORLANDO, FLORIDA 32804

PHONE: (800) 425-0012
FAX: (407) 425-5554 E-MAIL: bolt@netpass.com



The BALDWIN Boot

The *Baldwin* brace has been designed with a variety of Ankle and Foot pathologies in mind. Several of the conditions which are addressed with the *Baldwin* include painful arthritis of the ankle joint, acquired pes plano valgus (adult flat foot), posterior tibialis, and deformity secondary to Charcot and neuroarthropathy. The *Baldwin* Boot brace has proven most successful in the post-operative management of foot and ankle trauma. As an adjunct to open reduction and internal fixation, the *Baldwin* is second to none in its ability to offer edema control and stability of the affected joint, as well as, soft tissue protection.

The *Baldwin Boot* brace is fabricated with polypropylene and leather from a mold of the patient's limb. Casting techniques are crucial for successful application of the Baldwin Boot brace and are determined by the pathology presented.

For more information about the *Baldwin Boot* brace and how it can be used for your patients contact **BOLT Systems, Inc.** at 1-800-425-0012.